




G T.K Medical

WWW.TK-MEDICAL.COM



Patent Pending: PCT/CN2004/001174

# TKBAG®

The TKBAG/Retriever Bags are used in all types of laparoscopic surgery for removing excised tissues/specimen through a small surgical incision. They can be introduced through 5mm or 10mm or larger size trocar cannula. They comes in a variety of sizes (as shown in the table below) and is suitable for a wide range of laparoscopic procedures.

## Safety

Extra-durable bag  
prevents leaks and tears.  
Soft ring design reduces  
tissue damage risks.

## Quality

Strict, extensive raw material  
selection and 31 quality inspection  
steps on every TKBAG.

## Convenience: Easy 4 U.


Easy in, easy out, easy open, easy close; **Saves time!**  
Multiple sizes including extra-large.

## SPECIMEN RETRIEVAL BAG


Model	Product Code	Diameter (mm)	Bag Size (mm)	Volume (ml)
Standard	HSD-40	10	40	40
	HSD-60	10	60	75
	HSD-80	10	80	200
	HSD-100	10	100	350
	HSD-130	10	130	700
A	HSD-40/A	/	40	40
	HSD-60/A	/	60	75
	HSD-80/A	/	80	200
	HSD-100/A	/	100	350
	HSD-130/A	/	130	700
C	HSD-40/C	10	40	40
	HSD-60/C	10	60	75
	HSD-80/C	10	80	200
	HSD-100/C	10	100	350
	HSD-130/C	10	130	700
	HSD-130L/C	10	130	1000
S	HSD-40/S	5	40	40
	HSD-60/S	5	60	60



Model Standard



Model C



Model S



Model A



TÜVRheinland®

CE 0197

ISO13485

©2019 Guangzhou T.K Medical Instrument Co., Ltd.

Add: No. 7, Jianye Six Road, East Zone of Guangzhou Economic and Technological Development District, Huangpu District, 510765, Guangzhou, China

Mobile: +86135 3379 7435

Tel: +86 20 3761 6021/3229 0169

Fax: +86 20 37616020/ 32290395

Http: //www.tk-medical.com E-mail: lam@gtkmedical.com